



PRODUCT COMPARISON

performanceLine


UFSK
INTERNATIONAL

OSYS



Ophthalmic surgery



		<i>surgiTrend</i>		<i>surgiLine</i>		Status: 07/2026	
		Item No.: 2200 0012 / 2200 0010		Item No.: 2200 0013 / 2200 0008			
	Maximum dimensions (LxWxH) incl. armrests	655 x 620 x 1,150 mm 655 x 620 x 1,080 mm – (backrest and lifting column fully extended upwards)		Comfort seat Dynamic seat Formed seat	640 x 620 x 1,150 mm – 640 x 620 x 1,160 mm		
	Chassis	zinc undercoated, powder-coated steel frame					
	Chassis dimensions (LxW)	655 x 620 mm (incl. brake lever and foot switch)		640 x 620 mm (incl. brake lever and foot switch)			
	Castor	4x double castors 100 mm Ø 1x double castor 65 mm Ø		4x double castors 100 mm Ø			
	Brake system	central brake					
Battery	Battery box with lead gel batteries						
	Footswitch	height adjustment lifting column, with LED					
	Footswitch plate (LxWxH)	Item No.: 3351 1011 480 x 715 x 53 mm		Yes (optional)	Item No.: 3351 0014 480 x 715 x 53 mm		
	Traverse	zinc undercoated, powder-coated					
	Seating systems	Item No.: 2200 0012 Item No.: 2200 0010 –		Comfort seat Dynamic seat Formed seat	Item No.: 2200 0013 – Item No.: 2200 0008		
	Height adjustment sitting height	550 – 750 mm 600 – 800 mm –		Comfort seat Dynamic seat Formed seat	550 – 750 mm – 530 – 730 mm		
	Lifting height	200 mm, electromotive adjustment					
	Seat adjustment angle	-5°/+5° – –		Comfort seat Dynamic seat Formed seat	-5°/+5° – -5°/+5°		
	Adjustment angle backrest	-18°/+5° 29° –		Comfort seat Dynamic seat Formed seat	-18°/+5° – -18°/+5°		
	Autodynamic angle	seat – 12° –	backrest – 16° –	Comfort seat Dynamic seat Formed seat	seat – – –	backrest – – –	
	Upholstery (thickness)*	seat 65 – 95 mm 100 – 120 mm –	backrest 45 – 80 mm 95 mm –	Comfort seat Dynamic seat Formed seat	seat 65 – 95 mm – 35 – 60 mm	backrest 45 – 80 mm – 30 – 35 mm	
		Color card "skai Tundra" / "Stamskin Top"					
	Seat slide	No					

Our recommendations

Comfort seat: for temporal access to the operating field
Dynamic seat: for temporal access to the operating field
Formed seat: for superior approach to the operating field

* For upholstery, dimensional tolerances of ±5 mm must be taken into account.



PRODUCT COMPARISON

performanceLine

UFSK
INTERNATIONAL

OSYS



Ophthalmic surgery




surgiTrend
Item No.: 2200 0012 / 2200 0010

surgiLine
Item No.: 2200 0013 / 2200 0008

Status:
07/2026



Armrests	stable, ergonomic shape, supporting area (350 x 120 mm) height-adjustable, ball joint for full 360-degree range movement	
Height adjustment armrests	manual, asynchronous, non-sterile, stroke 320 mm	
Accessory armrest pair for neurosurgery	No	
Armrest pair with linear guide	Yes (optional): with one-hand release, supporting area 455 x 145 mm, adjustment range 110 mm, fully movable thanks to ball joint Item No.: 3351 1006	
Armrest pair with one-hand release	Yes (optional): with one-hand release, supporting area 455 x 145 mm, fully movable thanks to ball joint Item No.: 3351 1004	
Max. load	150 kg	



Our recommendations

Both surgeon's chairs are suitable for a team with only one surgeon.
If there is more than one surgeon, we recommend the one-handed armrest (Item No.: 3351 1004) and linear guide (Item No.: 3351 1006) for quick individual adjustment of the armrests.